



# OSTRACODTOXKIT F

## Test procedure



# 1

## **PREPARATION OF STANDARD FRESHWATER**

- VOLUMETRIC FLASK (1liter)
- VIALS WITH SOLUTIONS OF  
CONCENTRATED SALTS
- DISTILLED (or deionized) WATER



2

POUR THE 5 VIALS  
WITH CONCENTRATED SALT SOLUTIONS  
IN  $\pm$  800 mL DISTILLED WATER,  
IN THE 1 LITER VOLUMETRIC FLASK



**3**

- FILL THE FLASK TO THE 1 LITER MARK
- AERATE FOR AT LEAST 15 MINUTES



# 4

## HATCHING OF OSTRACOD CYSTS

POUR THE CONTENTS OF ONE VIAL  
WITH CYSTS IN THE PETRI DISH



**5**

TO ENSURE THE TRANSFER OF ALL THE CYSTS, THE VIAL SHOULD BE RINSED TWICE WITH 1 mL STANDARD FRESHWATER



# 6

## **INCUBATION OF THE CYSTS**

INCUBATE THE PETRI DISH  
FOR 48 HOURS AT 25 °C  
UNDER CONTINOUS ILLUMINATION  
OF MIN. 3 000 – 4 000 LUX



**7**

**PRE-FEEDING  
OF THE TEST ORGANISMS**

TAKE ONE VIAL  
WITH SPIRULINA POWDER  
AND FILL IT  
WITH STANDARD FRESHWATER



# 8




- SHAKE THE VIAL WITH THE SPIRULINA SUSPENSION
- POUR THE CONTENTS IN THE PETRI DISH CONTAINING THE HATCHED OSTRACODS AND SWIRL THE PETRI DISH GENTLY
- ALLOW THE OSTRACODS TO PRE- FEED FOR 4 HOURS



# 9

## **LENGTH MEASUREMENT OF FRESHLY HATCHED OSTRACODS**

PICK UP 10 OSTRACODS FROM THE  
HATCHING PETRI DISH  
WITH A GLASS MICROPIPETTE




**10**

TRANSFER THE OSTRACODS  
INTO ONE CUP  
OF THE MULTIWELL  
FOR "LENGTH MEASUREMENTS"



# 11


ADD ONE DROP OF LUGOL FIXATIVE  
TO THE CUP CONTAINING THE OSTRACODS  
AND WAIT UNTIL THE ORGANISMS ARE  
COMPLETELY IMMOBILE



**12**

PUT THE MICROMETER SLIP ON THE GLASS STAGE OF THE  
DISSECTION MICROSCOPE, IN THE CENTRE OF THE VISUAL FIELD

# 13



PUT THE MULTIWELL  
FOR LENGTH MEASUREMENTS  
ON THE STAGE OF THE DISSECTION  
MICROSCOPE, AND MEASURE  
THE LENGTH OF THE OSTRACODS

N.B : the smallest subdivisions of the  
micrometer slip are  $50\ \mu\text{m}$

FRESHLY HATCHED OSTRACODS HAVE  
A LENGTH OF ABOUT  $200\ \mu\text{m}$

# 14

SCORE THE LENGTH  
OF THE OSTRACODS  
ON THE “RESULTS SHEET”  
IN COLUMN “Day 0”





# 15

## **PREPARATION OF THE ALGAL FOOD SUSPENSION**

TAKE ONE TUBE WITH ALGAL BEADS  
AND POUR OUT THE STORAGE MEDIUM  
TAKING CARE NOT TO LOSE ANY BEAD  
DURING THE OPERATION





**16**

ADD 7 mL MATRIX DISSOLVING MEDIUM TO THE TUBE WITH ALGAL BEADS  
AND CLOSE THE TUBE WITH THE CAP



17

### DE-IMMOBILISATION OF THE ALGAE

SHAKE THE TUBE ON A VORTEX UNTIL THE MATRIX IN WHICH THE ALGAE ARE IMMOBILISED IS TOTALLY DISSOLVED AND THE ALGAE ARE SET FREE



CENTRIFUGE THE TUBE FOR 10 MINUTES  
AT 3000 rpm IN A CONVENTIONAL  
LAB CENTRIFUGE



CAREFULLY POUR OUT THE  
SUPERNATANT FROM THE TUBE

18



ADD 10 mL DISTILLED WATER  
TO THE TUBE WITH THE ALGAL PELLETT

**19**



CAP AND SHAKE THE TUBE  
TO RESUSPEND THE ALGAE



**20**

CENTRIFUGE THE TUBE AGAIN AT 3000 rpm FOR 10 MINUTES  
AND POUR OUT THE SUPERNATANT



ADD 10 mL STANDARD FRESHWATER  
TO THE TUBE WITH THE ALGAL PELLETT

**21**



CAP AND SHAKE THE TUBE  
TO RESUSPEND THE ALGAE



**22**

- POUR THE CONCENTRATED ALGAL SUSPENSION INTO A 25 mL VOLUMETRIC FLASK
- ADD STANDARD FRESHWATER TO THE 25 mL MARK
- CAP THE FLASK AND SHAKE TO OBTAIN A HOMOGENOUS ALGAL SUSPENSION




**23**

**ADDITION OF SEDIMENT, ALGAL FOOD  
AND OSTRACODS TO THE TEST PLATES**

PUT 2 mL STANDARD FRESHWATER  
INTO EACH WELL OF TWO TEST  
PLATES (multiwell for reference sediment  
and multiwell for test sediment)






# 24

## TEST PLATE FOR REFERENCE SEDIMENT

- FILL THE SPOON WITH REFERENCE SEDIMENT AND STRIKE OFF THE EXCESS SEDIMENT WITH THE PLASTIC SPATULA (the filled spoon then contains 500  $\mu$ l sediment)
- PUT 2 SPOONS (= 1000  $\mu$ l) REFERENCE SEDIMENT INTO EACH WELL OF THE TEST PLATE



# 25

## TEST PLATE WITH TEST SEDIMENT

- FILL THE SPOON WITH TEST SEDIMENT AND STRIKE OFF THE EXCESS SEDIMENT WITH THE PLASTIC SPATULA (the filled spoon then contains 500  $\mu\text{l}$  sediment)
- PUT 2 SPOONS (= 1000  $\mu\text{l}$ ) OF TEST SEDIMENT INTO EACH WELL OF THE TEST PLATE (use the tip of the spatula to perform the transfer)



## 26

- POUR THE ALGAL FOOD SUSPENSION FROM THE 25 mL FLASK INTO A BEAKER
- SHAKE THE BEAKER TO DISTRIBUTE THE ALGAE EVENLY
- PIPET 2 mL ALGAL SUSPENSION INTO EACH WELL OF THE TWO TEST PLATES



**27**

FILL THE LID OF THE HATCHING  
PETRI DISH WITH 10 mL STANDARD  
FRESHWATER



**28**

TRANSFER WITH THE GLASS  
MICROPIPETTE A NUMBER OF  
OSTRACOD NEONATES FROM THE  
HATCHING PETRI DISH INTO THE LID  
OF THIS DISH



**29**

TRANSFER 10 OSTRACODS FROM THE PETRI DISH LID INTO EACH WELL  
OF THE TWO TEST PLATES



**30**

### **INCUBATION OF THE TEST PLATES**

- COVER THE TWO TEST PLATES WITH A SHEET OF PARAFILM AND PUT THE LID ON TOP
- INCUBATE THE TEST PLATES AT 25 °C, IN DARKNESS, FOR 6 DAYS



**31**

**TEST SCORING – 1. TRANSFER OF THE OSTRACODS INTO A PETRI DISH**

WITH THE AID OF THE “LARGE MOUTH” MICROPIPET, SUCK UP PART OF THE SEDIMENT SUSPENSION FROM ONE CUP OF THE TEST PLATE AND TRANSFER IT INTO THE MICROSIEVE





## 32

- GENTLY RINSE THE CONTENTS OF THE MICROSIEVE WITH TAPWATER TO WASH OUT THE FINE SEDIMENT
- PROCEED FURTHER WITH THE STEPWISE TRANSFER OF THE SEDIMENT SUSPENSION TO THE MICROSIEVE AND RINSE EACH TIME THE CONTENTS OF THE MICROSIEVE



## 33

- ADD A FEW mL STANDARD FRESHWATER TO THE CUP
- MIX THE WATER WITH THE REMAINING SEDIMENT AND TRANSFER IT TO THE MICROSIEVE FOR RINSING.
- REPEAT THIS OPERATION UNTIL ALL THE SEDIMENT AND OSTRACODS HAVE BEEN TRANSFERRED INTO THE MICROSIEVE



## 34

- TURN THE MICROSIEVE UPSIDE DOWN AND WASH THE CONTENTS INTO A PETRI DISH, WITH THE AID OF A WASH BOTTLE CONTAINING TAPWATER
- REPEAT THE SEDIMENT TRANSFER AND RINSING OPERATIONS FOR ALL THE CUPS OF THE TWO TEST PLATES



**35**

**TEST SCORING – 2. MORTALITY SCORING**

PICK UP ALL THE LIVE OSTRACODS FROM THE PETRI DISH WITH A GLASS MICROPIPETTE  
AND TRANSFER THEM INTO ONE CUP OF THE “LENGTH MEASUREMENTS” MULTIWELL



# 36

- COUNT THE NUMBER OF LIVE OSTRACODS IN THE CUP
- SUBSTRACT THIS NUMBER FROM 10 (i.e. from the original number of ostracods put in the cup)
- SCORE THE OUTCOME (i.e. the number of dead ostracods) ON THE RESULTS SHEET
- REPEAT THIS OPERATION FOR ALL THE CUPS OF THE TWO TEST PLATES
- CALCULATE AND SCORE THE MEAN % OSTRACOD MORTALITY FOR ALL THE CUPS

# 37


## TEST SCORING – 3. LENGTH MEASUREMENT

*NB : only to be performed if*

*the percentage mortality is < 30%*

- ADD ONE DROP OF LUGOL FIXATIVE TO EACH CUP OF THE LENGTH MEASUREMENTS MULTIWELL WHICH CONTAIN THE LIVE OSTRACODS FROM THE TWO TEST PLATES





**38**

- WAIT UNTIL THE OSTRACODS ARE IMMOBILE
- MEASURE THE LENGTH OF EACH OSTRACOD FOLLOWING THE PROCEDURE INDICATED IN STEPS 13 & 14
- SCORE THE RESULTS IN THE CORRESPONDING “LENGTH” BOXES OF THE RESULTS SHEET

# OSTRACODTOXKIT F

## RESULTS SHEET

Mortality (D) and length ( $\mu$ ) of test organisms

Name of operator : TAMINGOL RIK

Date of performance of test : 16/07/2000

Test sediment : SAMPLE 304A

| LENGTH  | DAY 0 | REFERENCE SEDIMENT |      |      |      |      |      | TEST SEDIMENT |     |     |     |     |     |
|---|-------|--------------------|------|------|------|------|------|---------------|-----|-----|-----|-----|-----|
|   |       | DAY 6              |      |      |      |      |      | DAY 0         |     |     |     |     |     |
|   |       | Replicate          |      |      |      |      |      | Replicate     |     |     |     |     |     |
| Test organism                                 | 1     | 2                  | 3    | 4    | 5    | 6    | 1    | 2             | 3   | 4   | 5   | 6   |     |
| 1   | 200   | 1100               | 1100 | 1100 | 1100 | 1100 | 1100 | 160           | 150 | 170 | 160 | 170 | 150 |
| 2   | 200   | 1100               | 1100 | 1100 | 1100 | 1100 | 1100 | 150           | 150 | 160 | 160 | 160 | 150 |
| 3   | 200   | 1100               | 1000 | 1000 | 1000 | 1000 | 1100 | 150           | 160 | 160 | 160 | 160 | 160 |
| 4   | 200   | 1100               | 1100 | 1100 | 1100 | 1000 | 1000 | 150           | 160 | 160 | 160 | 160 | 160 |
| 5   | 200   | 1000               | 1100 | 1000 | 950  | 1000 | 1000 | 150           | 160 | 150 | 150 | 160 | 160 |
| 6   | 200   | 1000               | 1100 | 1100 | 950  | 950  | 900  | 150           | 160 | 150 | 160 | 160 | 160 |
| T   | 200   | 1000               | 1100 | 1100 | 1100 | 1100 | 1100 | 150           | 160 | 150 | 160 | 160 | 160 |
| II  | 200   | 300                | 300  | 300  | 350  | 350  | 300  | 500           | 600 | 600 | 600 | 650 | 550 |
| III   | 200   | 300                | 300  | 300  | 350  | 300  | 350  | 500           | M   | 600 | 600 | 650 | M   |
| 10  | 200   | 300                | 350  | M    | 350  | 300  | M    | 500           | M   | 550 | 550 | 600 | M   |
| <b>MORTALITY</b>                              |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Number of dead ostracods per replicate        |       | 0                  | 0    | 1    | 0    | 0    | 1    | 0             | 1   | 0   | 0   | 0   | 1   |
| Mean % mortality per replicate                |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean % mortality for all replicates           |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Standard deviation of mean % mortality        |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Variation coefficient of mean % mortality     |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| <b>GROWTH INHIBITION</b>                      |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean length of ostracods per replicate        |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean length for all replicates                |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean length increment per replicate           |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean length increment for all replicates      |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean growth inhibition (%) for all replicates |       |                    |      |      |      |      |      |               |     |     |     |     |     |
| Mean % growth inhibition to test sediment     |       |                    |      |      |      |      |      |               |     |     |     |     |     |

# 39

- PERFORM THE DATA TREATMENT OF THE RESULTS WITH AN APPROPRIATE PROGRAMME